

# INTERFERin® transfection reagent

## Short protocol - siRNA Transfection

### DAY 0: Cell seeding

→ Seed cells in **V** mL of serum containing medium according to the table below

Quantities per well, dish or flask

Culture vessel	Number of cells*	V = volume of medium during transfection
96-well	2 500 – 7 500	0.2 mL
24-well	15 000 – 35 000	1 mL
12-well	30 000 – 70 000	2 mL
6-well / 35 mm	100 000 – 200 000	4 mL
100 mm / flask 75 cm <sup>2</sup>	750 000 – 1.25 x 10 <sup>6</sup>	15 mL

\*For suspension cells, please refer to the complete protocol.

### DAY 1: Transfection = 1 nM siRNA

→ Perform transfection **in the presence of serum**


→ Transfect cells at **30-50% confluency**



Watch the video «*siRNA transfection using INTERFERin®*» on YouTube!

<http://www.youtube.com/watch?v=Yk8W8Cn0zjw>

Dilute **X** pmoles of siRNA  
in **W** µL of medium without serum  
Vortex 10 s and spin down



Add **Y** µL of  
INTERFERin® reagent

Vortex 10 s, spin down and  
incubate 10 min at RT

During the incubation time,  
remove **Z** mL of serum  
containing medium

Add transfection mix  
to the cells in serum  
containing medium

Incubate 24 to 72 h

Quantities per well, dish or flask

Culture vessel	W = volume of medium without serum	X = amount of siRNA added (1 nM)	Y = volume of INTERFERin® reagent	Z = volume of serum containing medium
96-well	50 µL	0.17 pmoles (2.4 ng)	0.75 ± 0.5 µL	0.125 mL
24-well	100 µL	0.6 pmoles (8.4 ng)	2 ± 1 µL	0.5 mL
12-well	200 µL	1.2 pmoles (17 ng)	4 ± 2 µL	1 mL
6-well / 35 mm	200 µL	2.2 pmoles (31 ng)	8 ± 4 µL	2 mL
100 mm / flask 75 cm <sup>2</sup>	500 µL	10.5 pmoles (147 ng)	40 ± 10 µL	10 mL

### DAY 2-3: Analyze gene silencing

See back page for optimization tips

Download complete protocol on <http://www.polyplus-transfection.com/resources/product-literature/>

# INTERFERin<sup>®</sup> transfection reagent

## Short protocol - Optimization Tips

### + Protocol Optimization

- + The siRNA final concentration may range from 1 to 50 nM depending on the cells and the target gene.
- + Cell confluency: between 30 and 50% at the time of transfection.
- + Check our online Cell Transfection Database at:  
<http://www.polyplus-transfection.com/resources/cell-transfection-database/>



### + Tips to increase cell viability of sensitive cells

- + Replace medium 4 h after transfection.
- + Check that silencing the target gene does not affect cell viability.

### + Use appropriate controls

- + Positive control: housekeeping gene (GAPDH or HPRT).
- + Negative control: mismatch, scramble or non-targeting sequence.
- + Be cautious with fluorescently labeled siRNA: 20 to 30 nM are needed to detect a signal, while only 1 nM can be sufficient for efficient silencing using INTERFERin<sup>®</sup>.

### + Good siRNA Transfection Practices

- + Store appropriately INTERFERin<sup>®</sup> ( $5 \pm 3^\circ\text{C}$ ). Do NOT freeze INTERFERin<sup>®</sup>.
- + Ensure that cells have been passaged more than twice and less than 20 times prior to transfection.
- + Discard overconfluent cells.
- + Check the half-lives of the protein and mRNA of interest and measure gene silencing accordingly 24 to 96 h after transfection.
- + Regularly check for mycoplasma contaminations.
- + Serum quality may drastically affect transfection efficiency. When purchasing a new batch of serum or trypsin, check cell viability as well as transfection efficiency.

**Note :** INTERFERin<sup>®</sup> is recommended for siRNA transfection. Please refer to the complete protocol available online at: <http://www.polyplus-transfection.com/resources/product-literature/>. Please use jetPRIME<sup>®</sup> for DNA/siRNA cotransfection.