

jetMESSENGER[®] transfection reagent

Short protocol - mRNA Transfection

DAY 0: Cell seeding

→ Seed cells in **V** ml of serum containing medium according to the table below

Quantities per well, dish or flask

Culture vessel	Number of cells*	V = volume of growth medium for cell seeding
96-well	7 500 - 12 500	0.125 mL
24-well	12 000 - 50 000	0.5 mL
12-well	80 000 - 100 000	1 mL
6-well / 35 mm	150 000 - 200 000	2 mL
100 mm / flask 75 cm ²	2 x 10 ⁶ - 4 x 10 ⁶	10 mL

*For specific cell type or suspension cells, please refer to the complete protocol.

DAY 1: Transfection


→ Perform transfection in the **standard cell growth medium**

→ Use jetMESSENGER[®] mRNA buffer only

→ Transfect cells at **60-80% confluency**

Dilute **X** µg of mRNA
in **W** µL of mRNA buffer
Vortex 10 s and spin down

Add transfection mix to the cells
in serum containing medium



Quantities per well, dish or flask

Culture vessel	W = volume of mRNA buffer	X = amount of mRNA added	Y = volume of jetMESSENGER [®] reagent
96-well	12.5 µL	0.1 µg	0.25 µL
24-well	50 µL	0.5 µg	1 µL
12-well	100 µL	1 µg	2 µL
6-well / 35 mm	200 µL	2 µg	4 µL
100 mm / flask 75 cm ²	1000 µL	10 µg	20 µL

DAY 2-3: Measure gene expression

See back page for optimization tips

Download complete protocol on <http://www.polyplus-transfection.com/resources/product-literature/>

jetMESSENGER[®] transfection reagent

Short protocol - Optimization Tips

+ Protocol Optimization

- + Test different mRNA amounts between 0.5X and 2X.
- + Test different mRNA/jetMESSENGER[®] ratios, 1:2 to 1:3.
- + Check our online Cell Transfection Database for cell specific protocols at:
<http://www.polyplus-transfection.com/resources/cell-transfection-database/>



Quantities per well, dish or flask

Culture vessel	W volume of mRNA buffer	X amount of mRNA added	Y volume of jetMESSENGER [®] reagent
96-well	12.5 µL	0.1 ± 0.05 µg	0.25 ± 0.05 µL
24-well	50 µL	0.5 ± 0.1 µg	1 ± 0.2 µL
12-well	100 µL	1 ± 0.2 µg	2 ± 0.4 µL
6-well / 35 mm	200 µL	2 ± 0.5 µg	4 ± 0.8 µL
100 mm / flask 75 cm ²	1000 µL	10 ± 2.5 µg	20 ± 4 µL

+ Tips to increase cell viability of sensitive cells

- + Wash cells 4 h after transfection.
- + Ensure that the mRNA is diluted in the mRNA buffer provided by Polyplus-transfection[®].
- + Analyze transfection at an earlier time point (e.g. at 24 h instead of 48 h).
- + Decrease the amount of mRNA added per well.
- + Decrease the volume of jetMESSENGER[®] reagent.
- + Use more stable chemically modified mRNA.
- + Check if the expressed protein may cause toxicity. If this is the case, reduce the amount of mRNA.

+ Good mRNA Transfection Practices

- + Store appropriately jetMESSENGER[®] (5 ± 3°C) and the mRNA (- 80°C).
- + Ensure that the quality of your mRNA is optimal. Preferably use mRNA purchased from an oligo supplier, instead of homemade transcribed mRNA.
- + Use a common reporter gene-encoding mRNA as a positive control (ex: Luciferase or GFP).
- + Ensure the medium is permissive to the transfection.
- + The use of chemically modified mRNA (Pseudouridine, 5' Methylcytosine, 5-methoxyuridine, etc...) could improve the transfection efficiency.
- + Ensure that all reagents are RNase-free.

Note: For more information regarding experimental conditions, please refer to the complete protocol available online at: <http://www.polyplus-transfection.com/resources/product-literature/>