

# Optimized PEI-mediated production of clinical grade viral vectors, PEIpro® and PEIpro®-HQ




Pascale Belguise, Valerie Keding, Jelena Vjetovic, Alengo Nyamay'Antu, Patrick Erbacher  
Polyplus-transfection, Bioparc, 850 Boulevard S. Brant, 67400 Illkirch, France


## Abstract

Gene- and cell therapy-based medicines are experiencing resurgence due to the advances made in developing new viral vector systems guided by safety, specificity and potency considerations. Adeno-Associated Virus (AAV) and Lentiviruses are commonly used in therapeutics and often produced using PEI-mediated transient transfection in HEK-293 or HEK-293T cells. The critical raw materials needed for cGMP vector production must be sourced from approved suppliers and should have gone through a rigorous testing program to reduce the risk of introducing adventitious agents into the production process. PEIpro® and PEIpro®-HQ, provided and manufactured by Polyplus-transfection®, are respectively the unique PEIs developed for transfection and suitable for virus production from process development up to use in GMP processes.

## Optimized PEI for R&D to Clinical-Grade Virus Production



**Optimization process of PEI polymer chemistry.** Whereas long polymer fragments lead to cell toxicity and short fragments lead to lower complexation potential (in red), optimized PEI size with a low polydispersity index decreases toxicity, while increasing complexation potential (in blue) and reproducibility in transfection.



**Manufacturing process of PEIpro® and PEIpro®-HQ reagents.** The linear structure of PEIpro® and the manufacturing process developed by Polyplus-transfection® ensure a high, stable and reproducible amount of protonable amines available for transfection while providing a fully deacylated molecule and an extremely low polymer chain length variation.

	PEIpro®	PEIpro®-HQ	
<b>Characteristics</b>			
1 mg/ml of Linear PEI in water	✓	✓	
Fully synthetic	✓	✓	
Manufacturing Process: Identical			
<b>Quality Controls</b>			
<b>Identity</b>	Molecular Weight of intermediate product	✓	✓
	Side chain content		✓
	Assay		✓
	Color, Clarity		✓
	pH		✓
<b>Purity</b>	Impurity profile		✓
	Endotoxin assay	✓	✓
<b>Safety</b>	Sterility test	✓	✓
	Mycoplasma	✓	✓
	Activity test	✓	✓
<b>Documentation</b>			
Certificate of Analysis	✓	✓	
Certificate of Origin (confirming the absence of components of animal origin)	✓	✓	
Detailed Batch Production Documentation to include in an IND or IMPD		✓	

**PEIpro®-HQ is a highly qualified grade of PEIpro® reagent.** The quality of PEIpro® and PEIpro®-HQ are continuously assessed during the manufacturing process with suitable control testing. In comparison to PEIpro®, a more extensive number of quality controls are performed on both the bulk material and the formulated product of PEIpro®-HQ to assess **identity, purity, safety, and potency.**


## Compatible with Various Production Culture Systems

Culture medium	Productivity using PEIpro®
DMEM (Gibco®)	+++
IMDM (Gibco®)	+++
Ham's F12 (Gibco®)	+++
BalanCD® HEK293 (Irvine Scientific®)	+++
FreeStyle™ 293 (Gibco®)	+++
FreeStyle™ F17 (Gibco®)	+++
HyClone™ HyCell™ TransFx™-H (GE Healthcare™)	+++
Pro293™ (Lonza®)	++
CD 293 (Gibco®)	-


**PEIpro® and PEIpro®-HQ are compatible with several serum-containing media and commercially available synthetic media for virus production in both adherent and suspension cells HEK-293 cells.**

FreeStyle™ is a trademark of Life Technologies™ Corporation.  
HyClone™, HyCell™ and TransFx™ are trademarks of General Electric Company.  
BalanCD® and Transfectory™ are trademarks of Irvine Scientifics  
PEIpro® is a trademark of Polyplus-transfection®.

## Reproducibility & Scalability




**Excellent lot-to-lot protein yield reproducibility using PEIpro®.** Suspension HEK-293 cells were seeded at  $1 \times 10^6$  cells/mL in serum-free medium and transfected with PEIpro® following the standard protocol. IgG3-Fc production was assayed 48 h after transfection using protein G affinity quantification (HPLC).



**DNA-PEIpro® complex size is identical, independently from the volume of transfection mix preparation.** Complexes were prepared with a DNA concentration of 0.01 mg per mL of complex volume at a DNA:Reagent ratio of 1:4, either in 10 mL or 1 L. The size of the complexes was then measured every ten minutes using the Zetasizer Nanometer ZS (Malvern Instrument, Malvern, UK).


## Gold standard for High Virus Production Yields



**Lentivirus and AAV production in HEK-293T and HEK-293 cells grown in suspension in BalanCD® HEK293 (Irvine Scientific®).** HEK-293T (lentivirus) and HEK-293 (AAV) cells were thawed directly into each medium and passaged every 3 to 4 days before going into a 2L benchtop bioreactor. Cells were seeded and cultured for 3 days before being transfected by PEIpro® (Polyplus). For transfection, four plasmids were used for lentivirus and three plasmids were used for AAV. Lentiviral and AAV titer were measured 48 and 72 hours post transfection (Data kindly provided by Généthon).

Cell culture system	Vector	Cells	Titer
CS10® / CF10®	AAV	Adherent HEK-293, HEK-293T	$10^{11}$ - $10^{13}$ VG / ml
Fixed-bed bioreactor (iCELLIS®)	AAV	Adherent HEK-293T	$10^{14}$ - $10^{16}$ VG / ml
Shaker Flask	AAV	Suspension HEK-293, HEK-293T	$10^9$ - $10^{10}$ VP / ml
Bioreactor	AAV	Suspension HEK-293, HEK-293T	$0.8$ - $1.5 \times 10^9$ - $10^{10}$ VG / ml
10 cm dish/75 cm <sup>2</sup>	Lentivirus	Adherent HEK-293, HEK-293T	$1$ - $2 \times 10^8$ TU / ml
HYPERflask®/HYPERstack®	Lentivirus	Adherent HEK-293, HEK-293T	$1$ - $2 \times 10^8$ TU / ml
Shaker Flask	Lentivirus	Suspension HEK-293F, HEK-293T	$2 \times 10^7$ - $10^{10}$ VP / ml
Bioreactor	Lentivirus	Suspension HEK-293, HEK-293T	$10^7$ IG / ml

**PEIpro® is the reagent of choice for virus production runs in most cell culture systems in both adherent and suspension cells** Irrespective of the cell culture-based system and production scale, PEIpro® and PEIpro®-HQ have led to efficient viral vector yields superior to  $10^7$  IG/mL and  $10^9$  VG/mL, respectively for lentiviruses and AAV.



**PEIpro® gives higher virus titers than Calcium Phosphate.** Lentiviruses were produced in adherent HEK-293 cells grown in serum-free culture medium, using 15 µg DNA and 30 µL PEIpro® per 75 cm<sup>2</sup> flask. Viral titers were determined by flow cytometry of supernatants 48 h after transfection.

## Conclusion

### Advantages of PEIpro® and its higher quality grade PEIpro® -HQ

- A PEI optimized for transfection, suitable for virus production (and for protein production).
- Synthetic animal free reagent manufactured according to a well-established process.
- Robust product with a great lot-to-lot reproducibility and a long shelf life.
- Ideal for process development up to large-scale therapeutic viral vector production.
- Highest quality PEI available with extra Quality Controls (identity, potency, safety and purity) and supplied with extensive documentation, **PEIpro®-HQ: Ideal for use in GMP processes.**